

# DERMATOLOGY WOODS LAMP



The BASCO make Woods Lamp Portable Mode comes with Twin Ultra Violet tubes to induce fluorescence in Dermatoses and an additional extra. White tube for Magnification. Both the Lights can be operated separately.

Examination of the skin or Hair with Woods Lamp is a time-trusted technique to aid in the clinical diagnosis of some cutaneous disorders. Wood Light examination should be performed in a dark room, and total darkness is also satisfactory.

1. Dermatoses in the intertriginous areas
2. Patches on the scalp with scaling or partial hair loss (especially with broken hairs)
3. Pigmentary alterations
4. Diagnosis of burrows in scabies
5. Blisters and punctate erosions on the dorsum of hands, forearms
6. Tinea versicolor

## DERMATOLOGY PORTABLE PUVA : UVA/UVB

Wave Tested :

2 to 4 Joules - 400 NM (NANOMETER)

Treatment for :

Vitiligo / Lucoderma

Lamp - UVB 4W

Black Light - FL4BL

Precaution : Use goggle / Protective goggle (without see-through) to protect eyes while using in face / lips. Joules will start reducing gradually while using. However after some time, bulbs can be replaced to get once again the correct wave length.



MODEL-PUVA-501

**CAUTION**  
Avoid using with one tube.  
Replace immediately with new tube to avoid overload to the transformer