



# ELECTRO CAUTERY

## DERMATOLOGY, E.N.T. & GENERAL

### BYFRICATOR

**BYFRICATOR B-202**



#### **SPECIAL FEATURES**

- Offers a choice of monopolar and bipolar current.
- Portable Model attached with Brief case for safe handling.
- Intensity variable accurately (50 - 100 WATTS OUTPUT)
- For minor coagulation and cutting, useful for E.N.T. Surgeons, Dermatologists, General Surgeons and cosmetologists.
- Ideal for instant coagulation, desiccation and fulguration.

#### **EFFICIENCY**

- This method is preferred by surgeons by removal of warts Acne Keloid, Corns,
- Cysts, Chonocrood, Tumors, Vaginal Cysts. Moles, Fobroma, Small tumours, Conjunctivities Epilation (Electrolysis) etc.,
- Excellent cosmetics result and obtained with little or no side effects.
- Abnormalities are removed bloodlessly.

#### **RESULT**

- Unit gives amazing results in the removal of the growth on the surface and within accessible orifices. After dehydration the growth drops off naturally within a few days or by giving local anaesthesia immediately.
- New skin appears gradually and takes the tint of surrounding skin.